

OSTEOCOIL

OSTEO  **COIL**
Nitinol Active Compression System

OsteoCoil™ Nitinol Active Compression System

5.5mm and 7.0mm
Surgical Technique Guide




TABLE OF CONTENTS

SYSTEM OVERVIEW	2
1 INDICATIONS FOR USE	2
2 CONTRAINDICATIONS	2
3 AO PRINCIPLES	2
4 IMPLANT DESIGN RATIONALE	3
SURGICAL TECHNIQUE	4
5 PREOPERATIVE PLANNING	4
6 BONE PREPARATION	5
7 FRACTURE REDUCTION	5
8 GUIDE WIRE PLACEMENT	6
9 DETERMINE IMPLANT SIZE	7
10 DRILL FOR IMPLANT	8
11 TAP FOR IMPLANT	9
12 INSERT IMPLANT	10
13 FINAL PLACEMENT	11
14 IMPLANT REMOVAL (OPTIONAL)	12
PRODUCT INFORMATION	13
15 INSTRUMENTS	13
16 STERILE IMPLANTS	15
17 INSTRUMENT TRAY LAYOUT	16
APPENDIX	17
18 REFERENCES/DISCLAIMERS	17



SYSTEM OVERVIEW

1 Indications for Use

The OsteoCoil™ Nitinol Compression System is indicated for use in bone reconstruction, osteotomy, arthrodesis, joint fusion, fracture repair, and fracture fixation of the small bones for wrist, hand, and foot.


2 Contraindications

Do not use the OsteoCoil Nitinol Compression System in cases of:


- inadequate bone quantity or bone quality, including osteoporotic or osteopenic bone
- foreign body sensitivity to implant material
- acute localized infections
- patients with limited blood supply
- patients who are unwilling or incapable of complying with post-operative care instructions
- patients with closed or inadequate medullary canals

3 AO Principles^{1,2}


The OsteoCoil Nitinol Compression System adheres to the AO Principles of Fracture Management such as:




Fracture reduction to restore anatomical relationships



Fracture fixation providing absolute or relative stability



Preservation of blood supply



Early and safe mobilization

¹Muller ME, M Allgower, R Schneider, H Willenegger. Manual of Internal Fixation. 3rd ed. Berlin Heidelberg New York: Springer. 1991


²Ruedi TP, RE Buckley, CG Moran. AO Principles of Fracture Management. 2nd ed. Stuttgart, New York: Thieme. 2007



IMPLANT FEATURES

4 Implant Design Rationale

The proprietary compressive technology allows the implant to extend as it is implanted. As the implant is inserted, the helical shaft feature expands, causing an opposing compression force between the implant head and threads. Typically, a compression implant provides an initial compression force which is lost as soon as there is bone length loss due to bone remodeling. The unique OsteoCoil implant results in a compression force applied across a bony fusion site even with loss of length throughout bone resorption.



Note: The OsteoCoil Nitinol Compression System can be used for a variety of indications. A subtalar arthrodesis is illustrated to demonstrate a universal operative technique.



SURGICAL TECHNIQUE

5 Preoperative Planning

Use AP & Lateral radiographs to evaluate the geometry of the target indication.



Available Implant Sizes

The following OsteoCoil implants are available sterile-packaged. Implants are universal and are not right and left specific.

5.5mm OsteoCoil Implant

Item#	Length	Item#	Length
OC55040	40mm	OC55058	58mm
OC55042	42mm	OC55060	60mm
OC55044	44mm	OC55065	65mm
OC55046	46mm	OC55070	70mm
OC55048	48mm	OC55075	75mm
OC55050	50mm	OC55080	80mm
OC55052	52mm	OC55085	85mm
OC55054	54mm	OC55090	90mm
OC55056	56mm		

7.0mm OsteoCoil Implant

Item#	Length	Item#	Length
OC70040	40mm	OC70090	90mm
OC70045	45mm	OC70095	95mm
OC70050	50mm	OC70100	100mm
OC70055	55mm	OC70105	105mm
OC70060	60mm	OC70110	110mm
OC70065	65mm	OC70115	115mm
OC70070	70mm	OC70120	120mm
OC70075	75mm	OC70125	125mm
OC70080	80mm	OC70130	130mm
OC70085	85mm		



6 Bone Preparation


Reduce the target operative site, exposing the bone if necessary. If needed, distract and prepare the area for fusion by removing cartilage until there is exposed subchondral bone. Leave the overall contours of the bones intact. Once all cartilage is removed, the subchondral bone may also be fenestrated to promote bone fusion. Place any graft material if desired.



7 Fracture Reduction

Reduce and clamp the fracture, fusion, or osteotomy with the bone reduction clamps. Consider reduction clamp placement so that it does not interfere with future implant insertion. Alternatively, temporary wires may be used to facilitate fracture reduction. The following wires are available in the instrument tray:

Item#	Size	Tip
For 7.0mm: OGW023	2.3mm	Blunt/Trocar
For 5.5mm: OGW016	1.6mm	Blunt/Trocar






8 Guide Wire Placement

Instruments:

- For 7.0mm:**
OGW023 SS Guide Wire 2.3mm Blunt/Trocar, Active Compression
- For 5.5mm:**
OGW016 SS Guide Wire 1.6mm Blunt/Trocar, Active Compression

Manually compress the joint or fracture before placing a wire. Keep bones in proper orientation and under compression throughout procedure. Additional Guide Wires may be used to fixate the fusion site.

TIP: Verify wire placement under fluoroscopy in at least two planes. Additional oblique views may be helpful to visualize the wire trajectory. Ensure there is sufficient Guide Wire length to allow for the threaded tip of the implant to cross the fusion site entirely.



9 Determine Implant Size

Instruments:

- ODG001 Depth Gauge, Active Compression, 5.5mm & 7.0mm

Place the Universal Cannulated Depth Gauge onto the wire and rest firmly against the bone. Read the length from the calibrated laser mark on the wire.


Note: Implant length will extend during insertion. The table below indicates the maximum extension for each implant diameter. Extension should be accounted for when selecting an implant.

Example: When using a 5.5mm diameter implant, expect up to 3.5mm of extension. As the implants come in 2 mm increments from 40-60mm and then 5mm increments, please size down appropriately such that the fully extended length of the implant matches the reading from the depth gauge.

Note: When using the countersink, please take into account the depth of countersinking in implant placement measurements. As discussed with respect to the implant size when a counter sink is not used, the total extension of the implant must be addressed.

Example: If the depth gauge reads 60mm, and the plan is to countersink 2mm, given the 4mm extension of a 7.0mm diameter implant, a 54mm implant should be utilized.

Implant Diameter	Maximum Extension
7.0mm	4.0mm
5.5mm	3.5mm






10 Drill for Implant

Instruments:

For 7.0mm:	
OTP070	7.0mm Tissue Protector, Active Compression
ODB070	Cannulated Drill Bit for 7.0mm, Stepped, Active Compression
For 5.5mm:	
OTP055	5.5mm Tissue Protector, Active Compression
ODB055	Cannulated Drill Bit for 5.5mm, Stepped, Active Compression

The Tissue Protector should be used when drilling for the OsteoCoil implant. Ensure that the Tissue Protector is firmly against bone and place the Cannulated Drill Bit over the Guide Wire.

Slowly advance the Drill Bit until the desired depth is achieved. Drill depth can be read from the Drill Bit calibration lines on the back of the Tissue Protector. Confirm appropriate depth under flouroscopy.



11 Tap for Implant

Instruments:


For 7.0mm:	
OCT070	7.0 Tap, Active Compression
SHBLK	Ratchet Handle, Small Hudson
For 5.5mm:	
OCT055	5.5 Tap, Active Compression
OOORG	Ratchet Handle, AO

Tapping should be performed prior to implant insertion to ensure proper implant extension.

Attach the Cannulated Tap to the Screwdriver Handle. Place the Cannulated Tap over the Guide Wire and slowly advance to the desired depth. Tapping depth can be monitored using the calibrated lines marked on the Cannulated Tap.

Note: Tapping should be performed by hand only, not under power. Tapping is strongly recommended for succesful implantation of the Osteocoil 7.0mm or 5.5mm Active Compression implant.

Tip: For soft or osteopenic bone, consider partially tapping, saving the last 1-2cm of native bone for thread purchase.





12 Insert Implant

Instruments:

For 7.0mm:

OCD070 T25 Cannulated Driver, Active Compression
SHBLK Ratchet Handle, Small Hudson

For 5.5mm:

OCD055 T15 Cannulated Driver, Active Compression
OSD055 T15 Solid Driver, Active Compression
OOORG Ratchet Handle, AO

Attach the Driver Shaft to the Ratchet Handle. Slowly advance the OsteoCoil implant into bone by turning in a clockwise direction until the implant head is fully seated against the exterior bone surface. Monitor progress under fluoroscopy. Once flush contact is achieved, perform additional implant turns (1.5 turns for 7.0mm, 1.75 turns for 5.5mm) to activate the compression feature. Activation of the compression feature can be confirmed by observing an increase in distance between the implant sleeve and implant threads from before to after activation.

Please note: T25 is used for 7.0mm implant, T15 is used for 5.5mm implant.


Precautions:

If the implant does not advance under torsion, immediately stop and remove the implant.

Make sure the implant is advancing during insertion by periodically monitoring under fluoroscopy. Twisting without advancing can cause implant failure.


Do not insert the OsteoCoil implant using power. Only insert using the Ratchet Driver Handle by hand.

It is recommended that upon implant seating into the tapped canal, the guidewire is removed to as to prevent guide wire seizing on implant upon removal.



13 Final Placement

Confirm final implant placement under fluoroscopy.





14 Implant Removal (optional)


Instruments:

For 7.0mm:	
OCD070	T25 Cannulated Driver, Active Compression
For 5.5mm:	
OCD055	T15 Cannulated Driver, Active Compression
OSD055	T15 Solid Driver, Active Compression
OOORG	Ratchet Handle, AO

Attach the desired length Driver Shaft to the Ratchet Handle. Slowly rotate the OsteoCoil implant in a counterclockwise direction until the implant is fully removed.

Confirm all implants have been successfully removed using fluoroscopy.

Please note: T25 is used for 7.0mm implant, T15 is used for 5.5mm implant.



PRODUCT INFORMATION

15 Instruments

5.5mm Instruments



T15 Cannulated Driver, Active Compression
OCD055




T15 Solid Driver, Active Compression
OSD055



SS Guide Wire 1.6mm Blunt/Trocar,
Active Compression
OGW016



5.5mm Cannulated Countersink,
Active Compression
OCC055



Ratchet Handle, AO
OOORG



Cannulated Drill Bit for 5.5mm, Stepped,
Active Compression
ODB055



5.5 Tap, Active Compression
OCT055



5.5mm Tissue Protector, Active Compression
OTP055



7.0mm Instruments




T25 Cannulated Driver, Active Compression
OCD070




SS Guide Wire 2.3mm Blunt/Trocar,
Active Compression
OGW023



7.0mm Cannulated Countersink,
Active Compression
OCC070



Ratchet Handle, Small Hudson
SHBLK



Cannulated Drill Bit for 7.0mm, Stepped,
Active Compression
ODB070




7.0 Tap, Active Compression
OCT070



7.0mm Tissue Protector, Active Compression
OTP070

Shared Instruments



Depth Gauge, Active Compression
ODG001



Wire Pusher, Active Compression
OWP557

Custom



Small Hudson Adaptor
SHADP



16 Sterile Implants





5.5mm OsteoCoil Implant	
Item #	Length (mm)
OC55040	40
OC55042	42
OC55044	44
OC55046	46
OC55048	48
OC55050	50
OC55052	52
OC55054	54
OC55056	56
OC55058	58
OC55060	60
OC55065	65
OC55070	70
OC55075	75
OC55080	80
OC55085	85
OC55090	90



7.0mm OsteoCoil Implant	
Item #	Length (mm)
OC70040	40
OC70045	45
OC70050	50
OC70055	55
OC70060	60
OC70065	65
OC70070	70
OC70075	75
OC70080	80
OC70085	85
OC70090	90
OC70095	95
OC70100	100
OC70105	105
OC70110	110
OC70115	115
OC70120	120
OC70125	125
OC70130	130



17 Instrument Tray Layout



APPENDIX

18 References/Disclaimers

Testing data on file at Conventus Flower Orthopedics.

This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Conventus Flower Orthopedics products.

The medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product instructions for use.

Postoperative management is patient-specific and dependent on the treating professional’s assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.

Please also refer to the OsteoCoil Nitinol Compression System reprocessing instructions for cleaning and sterilization information.

[illegible]

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There is no text or other content on the page.



Caution: Federal law restricts these devices to sale by or on the order of a licensed physician.

© Copyright 2025 | Flower Orthopedics Corporation d/b/a Conventus Flower Orthopedics.
OsteoCoil™ is a trademark of Conventus Flower Orthopedics.
All Rights Reserved. FOC 1114 Rev A 10/25

Conventus Flower Orthopedics
100 Witmer Road, Suite 280, Horsham, PA 19044
877.778.8587

ConventusFlower.com