


READY-FOR-SURGERY™ SOLUTIONS

FOREFOOT

FEATURES

- Bend zones facilitate optimal contouring
- Variable screw placement up to 20° cone
- Compression slots allow for fragment compression
- Low-profile, 1.5mm thickness
- Screws designed to sit flush within the plate



The Flower forefoot plate portfolio contains low-profile stabilization plates, designed to provide robust fixation while reducing implant prominence. Anatomic plate configurations with compression slots and all locking plates are available in multiple lengths and sizes to treat a wide variety of fractures and osteotomies.

ORDERING INFORMATION

Compression T-Plate, Straight

PRODUCT NO.	HOLES
FFP 317	7
FFP 318	8




Compression Y-Plates

PRODUCT NO.	HOLES
FFP 306	6
FFP 307	7



Compression T-Plate, Oblique

PRODUCT NO.	ORIENTATION	HOLES
FFP 107	Left	7
FFP 207	Right	7



Akin

PRODUCT NO.	HOLES
FFP 302	2



HOW DOES FLOWER ORTHOPEDICS BUILD A MORE EFFICIENT CASE?

FLOWERCUBES™

FlowerCubes™ are pre-packaged and ready for use, eliminating costly pre-op sterilization. The FlowerCube contains all of the implants and instruments required for specific surgical indications.



IMPLANT PORTFOLIO

Our full line of Foot & Ankle solutions covers forefoot to ankle with plating and screw solutions. The consistent screw head design maximizes efficiency.



READY-FOR-SURGERY WITH

EVERYTHING YOU NEED

E-KIT ADVANCED™ + ONE ADDITIONAL KIT

Yes, you can do a full Foot & Ankle case with two sterile-packed instrument kits!

The E-Kit Advanced streamlines the surgical steps by providing multifunctional instruments for both plates and screws.

When combined with an additional sterile-packaged kit, such as a drill bit kit, you will have all the tools required for a procedure, eliminating operating room back table clutter.

The robust, single-use instruments are new every time. FlowerCubes eliminate the traditional wear and tear of instrumentation and post-op decontamination.

