

## READY-FOR-SURGERY™ SOLUTIONS

# HEADLESS COMPRESSION SCREWS

### FEATURES

- Head Design – Tapered design to facilitate fragment compression.
- Compression – Pitch differential between the leading and trailing thread maximizes and maintains compression.
- Aggressive Screw Tip – Self-drilling and self-tapping flutes for smooth insertion even in dense cortical bone.
- Ready-For-Surgery™ – All Headless Compression Screws are individually sterile packaged to reduce the potential of infection.



## ORDERING INFORMATION

PRODUCT NO.	SIZE
HCS 106	D:2.4mm x L:6mm
HCS 108	D:2.4mm x L:8mm
HCS 110	D:2.4mm x L:10mm
HCS 112	D:2.4mm x L:12mm
HCS 114	D:2.4mm x L:14mm
HCS 116	D:2.4mm x L:16mm
HCS 118	D:2.4mm x L:18mm
HCS 120	D:2.4mm x L:20mm
HCS 122	D:2.4mm x L:22mm
HCS 124	D:2.4mm x L:24mm
HCS 126	D:2.4mm x L:26mm
HCS 128	D:2.4mm x L:28mm
HCS 132	D:2.4mm x L:32mm
HCS 136	D:2.4mm x L:36mm
HCS 140	D:2.4mm x L:40mm
HCS 210	D:3.0mm x L:10mm
HCS 212	D:3.0mm x L:12mm
HCS 214	D:3.0mm x L:14mm
HCS 216	D:3.0mm x L:16mm
HCS 218	D:3.0mm x L:18mm
HCS 220	D:3.0mm x L:20mm
HCS 222	D:3.0mm x L:22mm

PRODUCT NO.	SIZE
HCS 224	D:3.0mm x L:24mm
HCS 226	D:3.0mm x L:26mm
HCS 228	D:3.0mm x L:28mm
HCS 230	D:3.0mm x L:30mm
HCS 232	D:3.0mm x L:32mm
HCS 234	D:3.0mm x L:34mm
HCS 236	D:3.0mm x L:36mm
HCS 238	D:3.0mm x L:38mm
HCS 240	D:3.0mm x L:40mm
HCS 245	D:3.0mm x L:45mm
HCS 250	D:3.0mm x L:50mm
HCS 414	D:4.0mm x L:14mm
HCS 416	D:4.0mm x L:16mm
HCS 418	D:4.0mm x L:18mm
HCS 420	D:4.0mm x L:20mm
HCS 422	D:4.0mm x L:22mm
HCS 424	D:4.0mm x L:24mm
HCS 426	D:4.0mm x L:26mm
HCS 428	D:4.0mm x L:28mm
HCS 430	D:4.0mm x L:30mm
HCS 432	D:4.0mm x L:32mm
HCS 434	D:4.0mm x L:34mm

PRODUCT NO.	SIZE
HCS 436	D:4.0mm x L:36mm
HCS 438	D:4.0mm x L:38mm
HCS 440	D:4.0mm x L:40mm
HCS 445	D:4.0mm x L:45mm
HCS 450	D:4.0mm x L:50mm
HCS 455	D:4.0mm x L:55mm
HCS 460	D:4.0mm x L:60mm
HCS 740	D:6.5mm x L:40mm
HCS 745	D:6.5mm x L:45mm
HCS 750	D:6.5mm x L:50mm
HCS 755	D:6.5mm x L:55mm
HCS 760	D:6.5mm x L:60mm
HCS 765	D:6.5mm x L:65mm
HCS 770	D:6.5mm x L:70mm
HCS 775	D:6.5mm x L:75mm
HCS 780	D:6.5mm x L:80mm
HCS 785	D:6.5mm x L:85mm
HCS 790	D:6.5mm x L:90mm
HCS 791	D:6.5mm x L:95mm
HCS 792	D:6.5mm x L:100mm

# HOW DOES FLOWER ORTHOPEDICS BUILD A MORE EFFICIENT CASE?

## FLOWERCUBES™

FlowerCubes™ are pre-packaged and ready for use, eliminating costly pre-op sterilization. The FlowerCube contains all of the implants and instruments required for specific surgical indications.



## IMPLANT PORTFOLIO

Our full line of Foot & Ankle solutions covers forefoot to ankle with plating and screw solutions. The consistent screw head design maximizes efficiency.



## READY-FOR-SURGERY WITH

# EVERYTHING YOU NEED

E-KIT ADVANCED™ + ONE ADDITIONAL KIT

Yes, you can do a full Foot & Ankle case with two sterile-packed instrument kits!

The E-Kit Advanced streamlines the surgical steps by providing multifunctional instruments for both plates and screws.

When combined with an additional sterile-packaged kit, such as a drill bit kit, you will have all the tools required for a procedure, eliminating operating room back table clutter.

The robust, single-use instruments are new every time. FlowerCubes eliminate the traditional wear and tear of instrumentation and post-op decontamination.

