

READY-FOR-SURGERY™ SOLUTIONS

TIBIA PLATES

Flower tibia plates have a robust, anatomically-contoured design with ideal accommodation for the anatomy and pathology of complex tibia fractures.

The locking holes in the Flower tibia plates allow for variable angle placement within a 30° cone and accept both locking and non-locking screws.



ANTERIOR LATERAL DISTAL TIBIAL PLATE

- Holes: 7, 10
- Orientation: Left, Right
- Distal holes provide articular support and treat comminuted fractures



MEDIAL DISTAL TIBIAL PLATE

- Holes: 8, 10
- Orientation: Left, Right
- Pre-packaged with MIS insertion post for minimally invasive techniques limited incision techniques



DISTAL TIBIA A-PLATE

- Holes: 5
- Anterior or posterior placement



HOW DOES FLOWER ORTHOPEDICS BUILD A MORE EFFICIENT CASE?

FLOWERCUBES™

FlowerCubes™ are pre-packaged and ready for use, eliminating costly pre-op sterilization. The FlowerCube contains all of the implants and instruments required for specific surgical indications.



IMPLANT PORTFOLIO

Our full line of Foot & Ankle solutions covers forefoot to ankle with plating and screw solutions. The consistent screw head design maximizes efficiency.



READY-FOR-SURGERY WITH

EVERYTHING YOU NEED

E-KIT ADVANCED™ + ONE ADDITIONAL KIT

Yes, you can do a full Foot & Ankle case with two sterile-packed instrument kits!

The E-Kit Advanced streamlines the surgical steps by providing multifunctional instruments for both plates and screws.

When combined with an additional sterile-packaged kit, such as a drill bit kit, you will have all the tools required for a procedure, eliminating operating room back table clutter.

The robust, single-use instruments are new every time. FlowerCubes eliminate the traditional wear and tear of instrumentation and post-op decontamination.

