


## READY-FOR-SURGERY™ SOLUTIONS

# TIBIA PLATES

Flower tibia plates have a robust, anatomically-contoured design with ideal accommodation for the anatomy and pathology of complex tibia fractures.

The locking holes in the Flower tibia plates allow for variable angle placement within a 30° cone and accept both locking and non-locking screws.



### ANTERIOR LATERAL DISTAL TIBIAL PLATE

- Holes: 7, 10
- Orientation: Left, Right
- Distal holes provide articular support and treat comminuted fractures



### MEDIAL DISTAL TIBIAL PLATE

- Holes: 8, 10
- Orientation: Left, Right
- Pre-packaged with MIS insert on post for minimally invasive techniques limited incision techniques



### DISTAL TIBIA A-PLATE

- Holes: 5
- Anterior or posterior placement



