



A Revolutionary New Device

CONSTRUCTOR®

Knotless Adjustable
Button Technology



Unique Locking Design

The Constrictor® Rope Syndesmosis Repair Kit is indicated for the treatment of disruptions to the ankle syndesmosis ligament complex.

The technology provides surgeons the ability to treat syndesmotic injuries using a knotless adjustable button system comprising of UHMWPE suture and low profile titanium button.

BENEFITS

- Superior fixation under cyclic loading
- Physiologic stabilization of the ankle mortise



Distributed by Conventus Flower Orthopedics

DUNAMIS MEDICAL®


Make the switch to better design
dunamismedical.com

A Game Changer!


- Knotless repair
- Novel adjustable self-locking design to assist with graft tensioning

ADJUNCT TO FRACTURE REPAIR

- Reduce risk of rigid malreduction of ankle syndesmosis
- Button compatible with multiple plating systems
- Reduce rate of re-operations




Elongation During Cyclic Loading
(max load of 250N)



© 2021 Dunamis Medical. All Rights Reserved.



Load Required to Induce Failure - Cyclic Loading
(clinically defined as 3mm displacement or breakage)



© 2021 Dunamis Medical. All Rights Reserved.



Average maximum elongation during cyclic loading

Superior ultimate failure-to-load as compared to the market leader

*Test Data on file at Dunamis Medical

Constrictor[®] Rope Syndesmosis Repair Kit

Product No.

DLF 002

Constrictor Rope Syndesmosis Repair Kit - Contains 1 Constrictor Rope & 1 Drill Bit 3.7mm



Distributed by Conventus Flower Orthopedics

100 Witmer Road, Suite 280
Horsham, PA 19044

T: 877.778.8587
FlowerOrtho.com

FOC 1065 Rev B 12/21

DUNAMIS[®] MEDICAL[®]

Make the switch to better design
dunamismedical.com

© 2021 Dunamis Medical. All Rights Reserved.

