



READY-FOR-SURGERY™ FLOWERDERM™

ACELLULAR DERMAL ALLOGRAFTS

FlowerDerm™ is an acellular dermal allograft intended to be used as a skin substitute for wound covering. FlowerDerm™ is processed via proprietary techniques and is pre-hydrated which allows for improved biomechanical strength compared to other dermal allografts on the market. FlowerDerm™ Meshed form is intended to be used as a skin substitute for wound covering applications (eg. Diabetic Foot Ulcer, Venous Foot Ulcer).

READY-FOR-SURGERY™

FlowerDerm™ Acellular Dermal Allografts are pre-hydrated (no need for hydration prior to use), pliable and stored at room temperature (15°- 30°C)

SAFETY


Biocompatible, Sterility Assurance Level (SAL) of 10⁻⁶

EXTRACELLULAR MATRIX

Provides a scaffold for cellular attachment and proliferation as well as vascular in-growth

RAPID HEALING

Readily available supply of proteins (collagen, elastin, laminin & fibronectin) and growth factors to facilitate healing



ORDERING INFORMATION

PRODUCT NO.	SIZE	THICKNESS
FDM 555	5x5cm, Meshed	0.4mm – 1.0mm
FDM 558	5x8cm, Meshed	0.4mm – 1.0mm

Flower Orthopedics Corporation
100 Witmer Road, Suite 280
Horsham, PA 19044

T: 877.778.8587
F: 215.394.8904

FlowerOrtho.com