

# FlowerGraft™ Preshaped Bone Wedge System


FlowerGraft™ Preshaped Bone Wedge System is a complete line of preshaped structural allografts for Evans and Cotton procedures. FlowerGraft Preshaped Bone Wedges offer surgeons convenient, safe, reliable, Ready-for-Surgery™ options for performing reconstruction procedures in the foot and ankle.

<b>Surgically Efficient</b>	Allograft bone wedges are provided sterile, stored at ambient temperature, packaged in saline and do not require hydration prior to use.
	The Ready-for-Surgery™ procedure kits contain all the instruments necessary to select & insert the appropriate graft
	Anatomic preshaped wedges do not require intra-operative shaping.
<b>Smart Design</b>	Bicortical faces provide structural support and cancellous core provides an excellent scaffold for bone growth.
<b>Plating Options</b>	Low profile, procedure specific plates provide unique anatomic fit and accommodate the range of available graft options.


# FlowerGraft™ Preshaped Evans & Cotton Wedges

The Evans and Cotton Wedges are available in the following sizes to accommodate anatomic variability.

Part #	Description	H	H1	L	W
AGB 106	FlowerGraft™ Pre-Shaped Evans Wedge	6mm	4mm	22mm	22mm
AGB 108	FlowerGraft™ Pre-Shaped Evans Wedge	8mm	5mm	22mm	22mm
AGB 110	FlowerGraft™ Pre-Shaped Evans Wedge	10mm	5mm	22mm	22mm
AGB 112	FlowerGraft™ Pre-Shaped Evans Wedge	12mm	5mm	22mm	22mm



Part #	Description	H	L	W
AGB 305	FlowerGraft™ Pre-Shaped Cotton Wedge	5mm	22mm	16mm
AGB 306	FlowerGraft™ Pre-Shaped Cotton Wedge	6mm	22mm	16mm
AGB 307	FlowerGraft™ Pre-Shaped Cotton Wedge	7mm	22mm	16mm
AGB 308	FlowerGraft™ Pre-Shaped Cotton Wedge	8mm	22mm	16mm



# Evans & Cotton Osteotomy Procedure Kits

The Evans Osteotomy Procedure Kit (AGB 901) and Cotton Osteotomy Procedure Kit (AGB 902) are provided sterile and consist of

- ▶ Trial Sizer – For selecting appropriate sized graft
- ▶ Tamp – For inserting graft into osteotomy site
- ▶ K-Wires (1.6mm dia) – For provisional fixation

## Evans Osteotomy Procedure Kit



## Cotton Osteotomy Procedure Kit



## FlowerADC™ Kit (optionally available instrument)


The FlowerADC™ Kit maybe used to Distract/Compress and maintain stability in osteotomy procedures.



## Evans Plates (optional)

3-Hole plates are provided sterile and are available in multiple sizes to account for anatomic variability.

Part #	Description	W
AGB 116	Evans 3-Hole Plate, Right	16mm
AGB 120	Evans 3-Hole Plate, Right	20mm
AGB 216	Evans 3-Hole Plate, Left	16mm
AGB 220	Evans 3-Hole Plate, Left	20mm



©2017 Flower Orthopedics. All rights reserved. FlowerCube, FlowerCarriage, FlowerDerm, FlowerGraft, IndicationCube, IndicationCarriage, FlowerADC, FlowerCart, Ready-for-Surgery, FlowerPatch, FlowerFlo and FlowerPot are trademarks of Flower Orthopedics. Please refer to package insert for description, indications, contraindications, warnings, precautions, and other important information.



Flower Orthopedics Corporation • 100 Witmer Road, Suite 280 • Horsham, PA 19044  
T: 877.778.8587 • F: 215.394.8904 • [www.flowerortho.com](http://www.flowerortho.com)

FOC-FLWRGRFT 05/17