

## READY-FOR-SURGERY™ SOLUTIONS

# OPENING WEDGE PLATE

### ADVANTAGES OF THE FLOWER OPENING WEDGE PLATES

#### PRECISE CORRECTION

Seven wedges sizes available to accommodate all anatomy & pathology

#### ANATOMIC PLATE DESIGN

Different degrees of correction on native joint structures while maintaining metatarsal length.

#### LOW-PROFILE

Design minimizes soft tissue irritation, potential.

#### BILATERAL DESIGN

Plates fit equally well on the left & right metatarsals.



PRODUCT NO.	SIZE	DEGREES OF CORRECTION <sup>1</sup>
FLOW 020	2mm	6
FLOW 030	3mm	9
FLOW 035	3.5mm	10.5
FLOW 040	4mm	12
FLOW 045	4.5mm	13.5
FLOW 050	5mm	15
FLOW 060	6mm	18

<sup>1</sup>Easley, Mark E., et al. "Hallux Valgus: Proximal First Metatarsal Osteotomies." International Advances in Foot and Ankle Surgery, edited by Amol Saxena, 1st ed., Springer-Verlag London Ltd., 2012, p. 16.

