

READY-FOR-SURGERY™ SOLUTIONS

OPENING WEDGE PLATE

ADVANTAGES OF THE FLOWER OPENING WEDGE PLATES

PRECISE CORRECTION

Seven wedges sizes available to accommodate all anatomy & pathology

ANATOMIC PLATE DESIGN

Different degrees of correction native joint structures while maintaining metatarsal length.

LOW-PROFILE

Design minimizes soft tissue irritation, potential.

BILATERAL DESIGN

Plates fit equally well on the left & right metatarsals.



PRODUCT NO.	SIZE	DEGREES OF CORRECTION ¹
FOW 020	2mm	6
FOW 030	3mm	9
FOW 035	3.5mm	10.5
FOW 040	4mm	12
FOW 045	4.5mm	13.5
FOW 050	5mm	15
FOW 060	6mm	18

¹Easley, Mark E., et al. "Hallux Valgus: Proximal First Metatarsal Osteotomies." International Advances in Foot and Ankle Surgery, edited by Amol Saxena, 1st ed., Springer-Verlag London Ltd., 2012, p. 16.

HOW DOES FLOWER ORTHOPEDICS BUILD A MORE EFFICIENT CASE?

FLOWERCUBES™

FlowerCubes™ are pre-packaged and ready for use, eliminating costly pre-op sterilization. The FlowerCube contains all of the implants and instruments required for specific surgical indications.



IMPLANT PORTFOLIO

Our full line of Foot & Ankle solutions covers forefoot to ankle with plating and screw solutions. The consistent screw head design maximizes efficiency.



READY-FOR-SURGERY WITH

EVERYTHING YOU NEED

E-KIT ADVANCED™ + ONE ADDITIONAL KIT

Yes, you can do a full Foot & Ankle case with two sterile-packed instrument kits!

The E-Kit Advanced streamlines the surgical steps by providing multifunctional instruments for both plates and screws.

When combined with an additional sterile-packaged kit, such as a drill bit kit, you will have all the tools required for a procedure, eliminating operating room back table clutter.

The robust, single-use instruments are new every time. FlowerCubes eliminate the traditional wear and tear of instrumentation and post-op decontamination.

