

# Flex-Thread™ Distal Fibula Intramedullary Nail System

The Flex-Thread Distal Fibula Nail System is an innovative solution for the treatment and fixation of fibula fractures and osteotomies.

## SYSTEM FEATURES

### Innovative Flex-Thread Design

- Flexible titanium nail shaft
- Atraumatic blunt nose tip design
- Self-contours to native anatomy
- Screw-like insertion assists in maintaining reduction

### Proximal Nail Fixation

- Provides secure point contact fixation
- Multiple thread diameters
- Self-tapping reverse cutting flutes


### Syndesmotic Versatility

- Nail accepts flexible (rope) or rigid (screw) fixation

### Distal Nail Fixation

- ZERO Profile Interlocking Screws
  - Divergent 2.7mm distal screw fixation
  - Screw head engages fibula cortex for improved fixation





### Atraumatic Blunt-Tip Nose

- Guides nail across fracture site
- Minimizes cortical disruption

### Self-Cutting & Self-Tapping Flutes

- No tapping required
- Decreases procedure time

### Reverse-Cutting Flutes

- Precision cut design
- Facilitates nail removal


### Syndesmotic Versatility

- Accepts 3.5mm ZERO Profile cortex screws or flexible rope fixation

### Low Profile Nail Design

- Anatomic fit
- Minimizes soft-tissue disruption
- 2.7mm ZERO Profile Interlocking Cortex Screws

Nail Diameter		Nail Length	
3.5mm	130mm	180mm	
4.5mm	130mm	180mm	
5.5mm	130mm	180mm	



Conventus Flower Orthopedics  
 100 Witmer Road, Suite 280, Horsham, PA 19044  
 877.778.8587

[ConventusFlower.com](http://ConventusFlower.com)